

Bausch + Lomb Retina Catalogue

Together we create a better retina surgery experience



CATARACT



LASER



RETINA

BAUSCH + LOMB
See better. Live better.

BAUSCH+LOMB RETINA

Together we create a better retina surgery experience

Table of Contents

- 4 Vitrectomy Solutions for Any Practice
(Hardware and Accessories)
- 12 Illumination
- 14 Laser
- 18 Instrumentation
- 31 Vitreoretinal Tamponades
- 32 Anterior Capabilities
- 39 Product Reference



BL14304
Stellaris® PC Next Generation
Combined system with an integrated
laser with 1 Xenon + 1 Xenon Mercury
lamp assembly

BL14334
Stellaris® PC Next Generation
Combined system (without the
integrated laser) with 1 Xenon + 1
Xenon-Mercury lamp assembly

BL2295 Integrated Foot Pedal

BL4390 Foot Control Battery

BL4391 Foot Control Wall Charger

BL4393 Foot Control Charging Cradle

BL4394 Foot Control Power Cable



BL4363 Stellaris® PC Zero Level Bottle Hanger -
Dual Bottle



 = 1

BL3270 Fragmentation Handpiece

 = 1



BL3170 Stellaris Phaco Handpiece

 = 1




BL3320 Xenon Lamp Replacement

 = 1

BL3321 Mercury Lamp Replacement

 = 1



BL5275 20G Fragmentation Needle 15°
Single Use

 = 6



STERILE 

BL5277 23G Fragmentation Needle 30°
Single Use

 = 6



STERILE 

BL2272 Laser Upgrade Kit

 = 1


BL4364 Stellaris PC/Stellaris Mayo Tray Adapter

 = 1




Please contact your local representative to order a Mayo Tray.


BL5520WX
20G Combined Pack
Pack Includes:
A 20G Vitrectomy cutter
B Scleral plug package - 19 & 20
C 20G Wide angle light pipe
D 20G MVR knife



BL5523WVX
23G Combined Pack
Pack Includes:
A 23G Vitrectomy cutter
B 23G Wide angle light pipe
C 23G Valved ESA system
D 23G Scleral plugs



BL5525WVX
25G Combined Pack
Pack Includes:
A 25G Vitrectomy cutter
B 25G Wide angle light pipe
C 25G Valved ESA system
D 25G Scleral plugs



- All the above packs come with the items below.
- 1 Cassette with AFI tube assembly
 - 2 Second aspiration line
 - 3 Irrigation line
 - 4 Drapes for Mayo tray, screen and remote control
 - 5 Cassette stopper
 - 6 Test chamber
 - 7 Male luer coupler
 - 8 Irrigation MICS sleeves x 2
 - 9 Priming cup
 - 10 Needle wrench
 - 11 Pre-piercing spike assembly
 - 12 Second stopcock



* Combined Procedural X Packs are also available in Gravity.
* Combined Procedural X Packs not compatible with Stellaris Elite™

Combined standard procedural packs with AFI (Air Forced Infusion)*

BL5520W
20G Combined Pack
BL5523WV
23G Combined Pack
BL5525WV
25G Combined Pack



* Combined Standard Procedural Packs are also available in Gravity.


Basic Pack for procedural customization*

BL5500
AFI Combined, without cutter, light probe, valved ESA




* Combined Basic Procedural Packs are also available in Gravity.


BL5420WX
20G Posterior Pack
Pack Includes:
A 20G Vitrectomy cutter
B Scleral plug package - 19 & 20
C 20G Wide angle light pipe
D 20G MVR knife



BL5423WVX
23G Posterior Pack
Pack Includes:
A 23G Vitrectomy cutter
B 23G Wide angle light pipe
C 23G Valved ESA system
D 23G Scleral plugs




BL5425WVX
25G Posterior Pack
Pack Includes:
A 25G Vitrectomy cutter
B 25G Wide angle light pipe
C 25G Valved ESA system
D 25G Scleral plugs



All the above packs come with the items below.

- 1 Cassette with AFI tube assembly
- 2 Second aspiration line
- 3 Irrigation line
- 4 Drapes for Mayo tray, screen and remote control
- 5 Cassette stopper
- 6 Priming cup
- 7 Male luer coupler
- 8 Pre-piercing spike assembly
- 9 Second stopcock



Posterior standard procedural packs with AFI (Air Forced Infusion)*

BL5420W
20G Posterior Pack
BL5423WV
23G Posterior Pack
BL5425WV
25G Posterior Pack



* Posterior Standard Procedural Packs are also available in Gravity.

Basic Pack for procedural customization*

BL5400
AFI Posterior, without cutter, light probe, valved ESA



* Posterior Basic Procedural Packs with AFI are also available in Gravity.

* Combined Procedural X Packs are also available in Gravity.
* Combined Procedural X Packs not compatible with Stellaris Elite™

NEW!

Cutters for Stellaris PC*

BL5620 20G Vitrectomy Cutter



 = 6
STERILE 

BL5623S 23G Vitrectomy Cutter



 = 6
STERILE 

BL5625S 25G Vitrectomy Cutter



 = 6
STERILE 

BL5627S 27G Vitrectomy Cutter Bi-Blade® 



 = 6
STERILE 

BL3394 Handle Extension




 = 12
STERILE 




Valved Trocars

BL5281 23G ESA Kit-valved (kit of 3)



 = 6
STERILE 

BL5282 25G ESA Kit-valved (kit of 3)



 = 6
STERILE 

NEW!

BL5240 27G ESA Kit-valved (kit of 3)



 = 6
STERILE 

BL5283 23G Single Trocar-Valved



 = 6
STERILE 

BL5284 23G Replacement Valves



 = 6
STERILE 

BL5285 25G Single Trocar-Valved



 = 6
STERILE 

BL5286 25G Replacement Valves



 = 6
STERILE 

BL5241 27G Single Trocar-Valved



 = 6
STERILE 

Valved Trocars

58.20V Valved ONE-STEP™ Port Entry System
20G. Polyimide cannulas with built-in valve



 = 6x3
STERILE 

58.23STHV Valved Stealth Port Entry System 23G
58.25STHV Valved Stealth Port Entry System 25G
Stainless steel cannulas with removable valve cap



 = 6x3
STERILE 

Valved ONE-STEP™ Stainless Port Entry System

58.23SV Valved ONE-STEP™
Port Entry System, stainless 23G
58.25SV Valved ONE-STEP™
Port Entry System, stainless 25G
58.27SV Valved ONE-STEP™
Port Entry System, stainless 27G
Stainless steel cannulas with removable valve cap



 = 6x3
STERILE 

Standard Trocars

58.23STH Stealth Port Entry System 23G
58.25STH Stealth Port Entry System 25G
Stainless steel cannulas



 = 6x3
STERILE 

58.20 ONE-STEP™ Port Entry System 20G
58.23 ONE-STEP™ Port Entry System 23G
Polyimide cannulas



 = 6x3
STERILE 

ONE-STEP™ Stainless Port Entry System

58.23S ONE-STEP™
Stainless Port Entry System, 23G
58.25S ONE-STEP™
Stainless Port Entry System, 25G
58.27S ONE-STEP™
Stainless Port Entry System, 27G
Stainless steel cannulas



 = 6x3
STERILE 

* Cutters for Stellaris® PC not compatible with Elite™ Combo.

Focal

BL5720 20G
BL5723 23G
BL5725 25G

 = 6
STERILE 



Wide Angle

BL5820 20G
BL5823 23G
BL5825 25G

 = 6
STERILE 



NEW!

Mid-Field

BL5823M 23G

 = 6
STERILE 



BL5825M 25G

 = 6
STERILE 



BL5827M 27G

 = 6
STERILE 



NEW!

Diffuse Wide Angle

BL 5823D 23G

 = 6
STERILE 



BL 5825D 25G

 = 6
STERILE 



Adjustable Chandelier

BL5234 23G
BL5235 25G

 = 5
STERILE 



Light Sources / PHOTON

PHOTON-2X Dual Port Light Source

 = 1



Replacement Bulb

PLAMP PHOTON
P2-PLAMP PHOTON II

 = 1



Illumination Adapters

Synergetics® Adapters

57.20 for use with Alcon Accurus, Xenon
57.26 for use with DORC Associate
57.27 for use with DORC Brightstar
57.30 for use with Alcon Constellation
57.40 for use with Bausch + Lomb Stellaris® PC



Synergetics® Endoilluminators*

Widefield Diffusion

56.21P 20G
56.21.23P 23G
56.21.25PS 25G

 = 12
STERILE 



Corona Shielded Widefield

56.12.23P 23G
56.12.25P 25G

 = 12
STERILE 



Midfield

56.22P 20G
56.22.25P 25G
56.22.27P 27G

 = 12
STERILE 



Vivid

56.20.23P 23G
56.20.25P 25G

 = 12
STERILE 



Focal

56.02P 20G
56.02.23P 23G
56.02.25PS 25G
56.02.27P 27G

 = 12
STERILE 



Synergetics® Chandeliers*

Awh Vivid

56.54.25P 25G

 = 12
STERILE 



Oshima Vivid

56.54.27P 27G

 = 12
STERILE 



Awh

56.50.25P 25G

 = 12
STERILE 



Oshima Dual

56.50.29P 29G

 = 12
STERILE 



Synergetics® Illuminated Infusion Chandeliers*

Standard

56.30P 20G
56.30.23P 23G
56.30.25P 25G

 = 12
STERILE 



High-Flow

56.34.25P 25G

 = 12
STERILE 



Synergetics® Multifunctional Illumination*

Chow Membrane Pick

56.07.23P 23G

 = 12
STERILE 



High-Flow Irrigating Endoilluminator

56.14 19G

 = 12
STERILE 



Aspirating Endoilluminator

56.24.23P 23G

 = 12
STERILE 



* Adaptor is required with Stellaris® PC

Stellaris® PC

With fully integrated laser

Stellaris® PC now features fully integrated laser. Upgrade compatible with existing Stellaris® PC models.

BL2272 Laser Upgrade Kit

- Laser module
- Laser mounting brackets
- E-stop and cabling
- Key switch and cabling (2) keys
- Laser drawer replacement cover
- Laser power input module
- Anterior module with laser support
- Laser to power input module cabling
- Laser optical patch cord
- CAN and power cabling
- Upgrade software kit
- External new labels Goggles



BL2295 Integrated New Design Footpedal



Laser Filters



BL3236 Zeiss Fixed ESF



BL3240 Zeiss 2 Position Filter ESF



BL3245 Moeller Weidel Fixed ESF



BL3241 Leica 2 Position Filter ESF



BL3247 Leica 4-way Filter



BL3237 Leica Fixed ESF



BL3239 Laser Filter ESF Mount Topcon



BL3242 Adaptor for Fitting of 2 Position Filters



BL3246 LIO Replacement Bulb



BL3234 Laser On/Off Keys (1 set)



BL3235 Laser Safety Goggles



- BL3244 Interlock Bypass Connector Kit
- Microscope filter override key
 - Microscope filter connector
 - Laser room interlock connector



BL2298 Dual LIO



BL2296 Auxiliary Footpedal for LIO Use



Standard Laser Probe

- BL5290 20G STD. Laser PROBE Straight
BL5293 23G STD. Laser PROBE Straight
BL5295 25G STD. Laser PROBE Straight



Curved Standard Laser Probe

- BL5290C 20G STD. Laser PROBE Curved
BL5293C 23G STD. Laser PROBE Curved
BL5295C 25G STD. Laser PROBE Curved



Illuminated Laser Probe

- BL5290L 20G ILL. Laser PROBE Straight
BL5293L 23G ILL. Laser PROBE Straight
BL5295L 25G ILL. Laser PROBE Straight



Curved Illuminated Laser Probe

- BL5290LC 20G ILL. Laser PROBE Curved
BL5293LC 23G ILL. Laser PROBE Curved
BL5295LC 25G ILL. Laser PROBE Curved



Standard Aspirating Laser Probe

- BL5290ASPH 20G Hard Tip Asp. Laser PROBE
BL5293ASPH 23G Hard Tip Asp. Laser PROBE



Soft-tip Aspirating Laser Probe

- BL5290ASP 20G Soft Tip Asp. Laser PROBE S
BL5293ASP 23G Soft Tip Asp. Laser PROBE S



Directional™ Laser Probes

Directional™ II Laser Probe

55.27 20G
55.27.23 23G
55.27.25 25G
55.27.27 27G

 = 12
STERILE 




Directional™ Laser Probe*

55.26 20G
55.26.23 23G
55.26.25, 55.26.25S 25G

 = 12
STERILE 

NEW!

55.26.27 27G



55.26.231 Laser Probe Directional™ I 23G
Direct connect to Quantel Vitra and Supra

55.26.251 Laser Probe Directional™ I 25G
Direct connect to Quantel Vitra and Supra

Extendable Directional™ Laser Probe

55.26E 20G
55.26.23E 23G
55.26.25E 25G

 = 12
STERILE 



Inverted Directional™ Laser Probe*

55.36.23E 23G
55.36.25E 25G

 = 12
STERILE 



Straight Laser Probes

Straight Laser Probe*

55.21 20G
55.21.23 23G
55.21.25S 25G
55.21.27 27G

 = 12
STERILE 



Curved Laser Probes

Flexible Tapered Laser Probe*

55.25.23 23G
55.25.25 25G

 = 12
STERILE 



TruCurve Curved Laser Probe*, 37° Curve

55.23 20G
55.23.23 23G
55.23.25 25G
55.23.27 27G

 = 12
STERILE 



Curved Laser Probe*, 22° Curve

55.12 20G

 = 12
STERILE 



Illuminated Laser Probes*

Curved Fixed Extended Illuminated Laser Probe, 37° Curve

55.69P 20G
55.69.23P 23G
55.69.25P 25G

 = 12
STERILE 



Curved Maxillum™ Laser Probe, 37° Curve

55.29 20G

 = 12
STERILE 



Straight Fixed Extended Illuminated Laser Probe

55.62.23P 23G
55.62.25P 25G

 = 12
STERILE 



Flexible Tapered Illuminated Laser Probe

55.70.23P 23G
55.70.25P 25G

 = 12
STERILE 



Illuminated Directional™ Laser Probe

55.47P 20G
55.47.23P 23G
55.47.25P 25G

 = 12
STERILE 



Illuminated Directional™ II Laser Probe

55.48.23P 23G
55.48.25P 25G

 = 12
STERILE 



Straight Maxillum™ Illuminated Laser Probe

55.22 20G

 = 12
STERILE 



Laser Source to which
Synergetics® Laser Probes will Connect:
Iridex GL / GLx / TX
Iridex IQ 532 / IQ 577
Alcon Eyelite
Alcon Purepoint
Alcon Constellation
Bausch + Lomb Millenium®
Bausch + Lomb Stellaris®PC
Ellex Solitaire

* Adaptor is required with Stellaris® PC

*Direct Connect models available for
Quantel Vitra Quantel Supra and Lumenis

Handle Designs Storz® Ophthalmic Forceps

- Fully disposable
- Ergonomic Handle with actuation from two sides
- Non reflective instrument shaft
- Color coded handles for gauge identification
- 20G, 23G, 25G, 27G



Selected sterile single use instruments are available as component of a Per Procedure Tray.
The Per Procedure Tray is a complete kit of sterile single use instruments selected by you. It is a flexible solution designed to provide safety, consistency and convenience in your operating room.
Assemble your individual kit of sterile single use instruments with the new online Per Procedure Tray Configurator.
For more information visit www.storzeye.eu



Endgripping Forceps

65023S-12 23G
65025S-12 25G

STERILE = 12

NEW!

SUI13-12 27G



Asymmetric Peeling Forceps

SUI03-12 23G
SUI05-12 25G

STERILE = 12



ILM Forceps

SUI09-12 20G
SUI07-12 23G
SUI08-12 25G

STERILE = 12



Serrated Forceps

SUI10-12 20G
SUI11-12 23G
SUI12-12 25G

STERILE = 12

NEW!

SUI14-12 27G



Membrane Picks

SUK01 23G

STERILE = 8



SUK02 25G

STERILE = 8



Bipolar Erasers

E7928 Scleral Bipolar Eraser,
Lemosa Type Connector

STERILE = 12



20G

E7929 20G Bipolar Eraser, Tapered,
Lemosa Type Connector

STERILE = 12



23G

E7930 23G Bipolar Eraser, Tapered,
Lemosa Type Connector

STERILE = 12



E7931 23G Bipolar Eraser, Tapered,
2-pin Connector

STERILE = 12



25G

E7924 25G Bipolar Eraser, Tapered,
2-pin Connector

STERILE = 12



E7925 25G Bipolar Eraser, Blunt,
2-pin Connector

STERILE = 12



E7926 25G Bipolar Eraser, Tapered,
Lemosa Type Connector

STERILE = 12



E7927 25G Bipolar Eraser, Blunt,
Lemosa Type Connector

STERILE = 12



NEW!

27G

E7932 27G Bipolar Eraser, Tapered,
Lemosa Type Connector

STERILE = 12



E7933 27G Bipolar Eraser, Tapered,
2-pin Connector

STERILE = 12



CX9430 Reusable Bipolar Cord
with Lemosa Connector

= 1

CX9400 Reusable Bipolar Cord with
2-pin Connector

= 1

Aspiration Handpieces

SUJ03 23G Passive Aspiration Handpiece. 1.5mm soft tip cannula

STERILE

= 12

SUJ02 25G Passive Aspiration Handpiece. 1.5mm soft tip cannula

STERILE

= 12

SUJ06 23G Passive Aspiration Handpiece. Blunt needle

STERILE

= 12

SUJ07 25G Passive Aspiration Handpiece. Blunt needle

STERILE

= 12

SUJ08 20G Backflush/Extrusion Handpiece. Blunt tip. Grey color code ring

STERILE

= 6

SUJ09 20G Backflush/Extrusion Handpiece. Soft tip brush. Grey color code ring. Not compatible with valved trocar/cannula system

STERILE

= 6

SUJ10 20G Backflush/Extrusion Handpiece. Soft tip. Grey color code ring. Not compatible with valved trocar/cannula system

STERILE

= 6

SUJ11 23G Backflush/Extrusion Handpiece. Blunt tip. Green color code ring.

STERILE

= 6

SUJ12 23G Backflush/Extrusion Handpiece. Soft tip brush. Green color code ring. Not compatible with valved trocar/cannula system

STERILE

= 6

SUJ13 23G Backflush/Extrusion Handpiece. Soft tip. Green color code ring. Not compatible with valved trocar/cannula system

STERILE

= 6

SUJ14 25G Backflush/Extrusion Handpiece. Blunt tip. Blue color code ring.

STERILE

= 6

NEW!

SUJ17 27G Backflush/Extrusion Handpiece. Blunt tip. Magenta color code ring.

STERILE

= 6

MVR Knives

LaseredgePlus®

LP7519 19G MVR Blade A:1.37mm

STERILE

= 6

LP7520 20G MVR Blade A:1.2mm

Handle Designs
Pinnacle 360°™

Pinnacle 360° is our high quality line of fully disposable, hand-held retinal surgical instruments.

- Offers wide variety of disposable tip options
- Fully disposable
- 360° actuation to approach the retina from any angle
- Color coded handles for gauge identification
- 20G, 23G, 25G, 27G, Foreign Body



Pinnacle 360°

20G

23G

25G

27G

Foreign Body

ILM Forceps

Eckardt Forceps

39.08PIN 20G
39.08.23PIN 23G
39.08.25PIN 25G
39.08.27PIN 27G

STERILE  = 5



End-Grasping Forceps

09.12PIN 20G
09.12.23PIN 23G
09.12.25PIN 25G

STERILE  = 5



Fine Tip Eckardt Forceps

09.10.23PIN 23G
09.10.25PIN 25G

STERILE  = 5



Microserrated™ End-Grasping Forceps

39.12.23PIN 23G
39.12.25PIN 25G

STERILE  = 5



Tano Asymmetrical Micro Forceps

39.03.23PIN 23G
39.03.25PIN 25G
39.03.27PIN 27G

STERILE  = 5



Heavy Membrane Forceps

Pyramid Tip Forceps

11.04.23PIN 23G
11.04.25PIN 25G

STERILE  = 5



Serrated Jaw Forceps

11.06.23PIN 23G
11.06.25PIN 25G

STERILE  = 5



Dual Purpose Forceps

Apte Corset Forceps

39.07.23PIN 23G
39.07.25PIN 25G

STERILE  = 5



Microserrated™ Tewari Omnicep Forceps

39.15.23PIN 23G
39.15.25PIN 25G

STERILE  = 5

NEW!

39.15.27PIN 27G



Foreign Body

BB Removal Forceps

11.40PIN 18G

STERILE  = 5



Foreign Body Forceps

11.20PIN 19G

STERILE  = 5



Foreign Body Magnet

12.21PIN 18G

STERILE  = 5



Pick Forceps

Pick Forceps

09.05.23PIN 23G
09.05.25PIN 25G

STERILE  = 5



Scissors

Horizontal Scissors Curved

23.76PIN 20G
23.76.23PIN 23G
23.76.25PIN 25G

STERILE  = 5



Morris Vertical Scissors

10.24PIN 20G
10.24.23PIN 23G
10.24.25PIN 25G

STERILE  = 5



Myopic Forceps

Fine Tip Eckardt Forceps

09.10.23PIN-M 23G
09.10.25PIN-M 25G
• 37mm

STERILE  = 5



Microserrated™ Tewari Omnicep Forceps

39.15.23PIN-M 23G
39.15.25PIN-M 25G
• 37mm

STERILE  = 5



I-Pack™

18060 COMPONENTS

STERILE  = 10

- 1mL Tuberculin Syringe
- 1 Pair Synthetic Latex-Free Gloves
- 130G 1/2" Needle
- 127G 1/2" Needle
- 2 18G 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 13.5-4mm Scleral Marker

18061 COMPONENTS

STERILE  = 10

- 2 1mL Tuberculin Syringes
- 1 Pair Synthetic Latex-Free Gloves
- 130G 1/2" Needle
- 127G 1/2" Needle
- 2 18G 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 13.5-4mm Scleral Marker
- 1 Eye Drape (16" x 16")

18064 COMPONENTS

STERILE  = 10

- 2 1mL Tuberculin Syringes
- 1 Pair Synthetic Latex-Free Gloves
- 130G 1/2" Needle
- 127G 1/2" Needle
- 2 18G 1" Needles
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 13.5-4mm Scleral Marker
- 1 Eye Drape (16" x 16")
- 1 Disposable Forceps, 4"
- 4 4" x 4" Gauze Pads

18067 COMPONENTS

STERILE  = 10

- 130G 1/2" Needle
- 1 Barraquer Eyelid Speculum, Adult
- 13.5-4mm Scleral Marker

18069 COMPONENTS

STERILE  = 10

- 1 Pair Synthetic Latex-Free Gloves
- 130G 1/2" Needle
- 6 Cotton Swabs 6"
- 1 Barraquer Eyelid Speculum, Adult
- 13.5-4mm Scleral Marker



Disposable Lenses

Awh Sutureless Lens Ring

41.00
• Used to stabilize a vitrectomy lens
• 3 sizes/pack

STERILE  = 10



Flat

VFD10
• Field of view: 36°
• Image magnification: 1.00x

STERILE  = 10



Bi-Concave

VBCD10
• Field of view: 25°
• Image magnification: 0.80x

STERILE  = 10



Magnifying

VMD10
• Field of view: 30°
• Image magnification: 1.50x

STERILE  = 10



Widefield

VWFD10
• Field of view: 48°
• Image magnification: 0.50x

STERILE  = 10



30° Prism

V30PD10
• Field of view: 33° (offset 30°)
• Image magnification: 1.00x

STERILE  = 10



Flat, SSV

VFLATSSVD
• Field of view: 30°
• Image magnification: 0.92x

STERILE  = 10



Single Use Bio Lenses

Single Use 20D Lens

V20LCD
STERILE  = 10

Single Use 28D Lens

V28LCD
STERILE  = 10

VOLK Optical Lenses

VCRLVITSSV/SVCRLVITSSV Central Retinal SSV
Field of view: 73°/88°
Image magnification: 0.71x



VMQVIT/SVMQVIT MiniQuad®
Field of view: 106°/127°
Image magnification: 0.39x



VMQVITACS/SVMQVITACS MiniQuad® ACS®
Peripheral indirect vitreoretinal procedures
Field of view : 106°/127°
Image magnification: 0.48x



VMQXLVIT/ SVMQXLVIT MiniQuad® XL
Field of view : 112°/134°
Image magnification: 0.39x



VQFL/SVQFL Quadr Aspheric®
Wide field examination and laser treatment
Field of view: 120°/144°
Image magnification: 0.51x
Laser spot: 1.97x



V20LC/SV20LC Binocular Lens 20D Large
Field of view: 46°/60°
Image magnification: 3.13x
Laser spot: 0.32x



V28LC Binocular Lens 28D Large
Field of view: 53°/69°
Image magnification: 2.27x
Laser spot: 0.44x



V30LC Binocular Lens 30D Large
Field of view : 58°/75°
Image magnification: 2.15x
Laser spot: 0.47x



VSQUAD160 SuperQuad® 160
Extreme wide field examination and laser treatment
Field of view : 160°/165°
Image magnification : 0.50x
Laser spot: 2.00x



VSFNC Super Field® NC Wide Field,
Pan retinal examination
Field of view : 95°/116°
Image magnification : 0.76x
Laser spot: 1.30x



VSPXL SuperPupil® XL,
Small pupil pan retinal examination
Field of view: 103°/124°
Image magnification: 0.45x
Laser spot: 2.20x



VG4HM 4 Mirrors Lens Gonio
G-4 High Mag
Field of view: 4x64°
Image magnification: 1.50x
Laser spot: 0.67x



V78C Double Aspheric – 78D
Field of view: 81°/97°
Image magnification: 0.93x
Laser spot: 1.08x



V90C 90D
Field of view: 74°/89°
Image magnification: 0.76x
Laser spot: 1.32x



VS66
Field of view: 80°/96°
Image magnification: 1.0x
Laser spot: 1.0x



VDGTLWF Digital Wide Field
Field of view: 103°/124°
Image magnification: 0.72x
Laser spot: 1.39x



VHRC HR Centralis
Field of view: 74°/88°
Image magnification: 1.08x
Laser spot: 0.93x



VVITHAN-LG VitreoLens Handle
Used with MiniQuad® and Central Retinal



VVITHAN-MXQL VitreoLens Handle
Used with HRX, MiniQuad® XL and Super Macula



VOLK Optical Lenses

VAC Area Centralis
Field of view: 70°/84°
Image magnification: 1.06x
Laser spot: 0.94x



VG4HMHAN2 3 Mirror High Mag Gonio



VHRXVIT HRX Vit
Field of view: 130°/150°
Image magnification: 0.43x



VHRWF High Resolution Wide Field
Field of view: 160°/165°
Image magnification: 0.50x
Laser spot: 2.0x



VCRLVIT Central Retinal
Field of view: 73°/88°
Image magnification: 0.71x



VHRXVITSSV HRX Vit SSV
Field of view: 130°/150°
Image magnification: 0.43x



V3MIR Mirror Laser Lens
Field of view: 60°/66°/76°
Image magnification: 1.06x
Laser spot: 0.94x



VMQVITSSV MiniQuad SSV



VMQVITSSVACS MiniQuad ACS SSV VIT
Field of view: 106°/127°
Image magnification: 0.48x



VG3 3 Mirror Lens Gonio



VDVVIT DynaView 156
Field of view: 95°/127°
Image magnification: 0.39°



VSRs2 Suture Ring
Field of view:
Image magnification:
Laser spot:



VG4LNF 4 Mirror Lens Gonio
Field of view: 60°/66°/76°
Image magnification: 1.06x
Laser spot: 0.94x



VMQXLVITSSV MiniQuad XL SSV
Field of view: 112°/134°
Image magnification: 0.39x



VSCA Sterilization Tray Small



DDMS™ Diamond Dusted Membrane Scrapers

The Diamond-dusted, soft, silicone tip has been designed to identify and lift the edge of the membrane quickly and easily for your forceps. Available in several styles, including a retractable version, the patented DDMS has become a staple in operating rooms around the world.

DDMS benefits:

- Aids in the atraumatic removal of immature membranes and pigmented cells from the retina
- Assists atraumatic separation of hyaloids from neurosensory retina
- Especially beneficial for relieving traction from the rim of macular holes or retinal breaks in an atraumatic fashion

Tano

20.04 20G
20.04.23 23G
20.04.25S (Stiff) 25G

NEW!

20.04.27 27G



STERILE = 5

Stiff Tano

20.07 20G



STERILE = 5

Retractable Tano

20.16 20G
20.16.23 23G
20.16.25 25G



STERILE = 3

STERILE = 5

STERILE = 5

Adjustable Tano with Pick

20.18 20G



STERILE = 3

Disposable Infusion Cannulas

6mm Tip

30.32 20G

- Ideal aphakic length



STERILE = 5

6mm Tip

30.32.23 23G

30.32.25 25G

- For severe vitreous hemorrhaging
- For use with Synergetics® steel cannula sets (sold separately)



STERILE = 6

Bipolar

Bipolar Cautery

53.00.23 23G

53.00.25 25G

NEW!

53.00.27 27G



STERILE = 6

Iris Retractors

Disposable

40.02-2

40.02-6

- 5 retractors/set in case one is lost or damaged



STERILE = 2

STERILE = 6

Visco Dissector

Hubbard Visco Dissector

12.06.25 25G 28G Tip

NEW!

12.06.27 27G

- Injects viscoelastic for delaminating large membranes from the retina



STERILE = 6

VRL100 DK-Line® 5ml

STERILE = 1

VRL110 DK-Line® 7ml

STERILE = 1



VRL200 Okta-Line® 5ml

STERILE = 1



VRL500 Oxane® 1300, 10ml Syringe

STERILE = 1

VRL600 Oxane® 5700, 10ml Syringe

STERILE = 1



VRL700 Oxane® HD, 10ml Syringe

STERILE = 1



CX5710 Viscous Fluid Extraction Kits

STERILE = 6



CX5720 Viscous Fluid Injection Kits for Oxane range

STERILE = 6



30.90 Silicone Oil collector chamber

STERILE = 6



NEW!

BL7620 High Flow VFI Cannula 20G

STERILE = 6



BL7623 High Flow VFI Cannula 23G

STERILE = 6



BL7625 High Flow VFI Cannula 25G

STERILE = 6



BL5110

Stellaris® Premium Vacuum
Phacoemulsification Pack

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Auxillary drape

C

Screen drape

D

Remote control drape

E

Irrigation sleeves x 2

F


Plug

G

Test chamber

H

Needle wrench



BL5111

Stellaris® Basic Vacuum
Phacoemulsification Pack

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Irrigation sleeve

C

Plug

D

Test chamber

E

Needle wrench



BL5112

Stellaris® StableChamber™ Vacuum Pack

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Auxillary drape

C

Screen drape

D

Remote control drape

E

Irrigation sleeves x 2

F

Plug

G


Test chamber

H

Needle wrench

I

StableChamber



BL5113

Stellaris® StableChamber™ Phaco Pack
with MICS Straight Needle

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Auxillary drape

C

Screen drape

D

Remote control drape

E

Irrigation sleeves x 2

F

Plug

G

Test chamber

H


Needle wrench

I

StableChamber

J

MICS straight needle



BL5113A

Stellaris® StableChamber™ Phaco Pack
with MICS Angled Needle

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Auxillary drape

C

Screen drape

D

Remote control drape

E

Irrigation sleeves x 2

F

Plug

G

Test chamber

H


Needle wrench

I

StableChamber

J

MICS angled needle



BL5114

Stellaris® Premium Vacuum Phaco Pack
with MICS Straight Needle

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Auxillary drape

C

Screen drape

D

Remote control drape

E

Irrigation sleeves x 2

F

Plug

G


Test chamber

H

Needle wrench

I

MICS straight needle



BL5114A

Stellaris® Premium Vacuum Phaco Pack
with MICS Angled Needle

=6

STERILE

Pack Includes:

A

VFM cassette assembly
and I/A set

B

Auxillary drape

C

Screen drape

D

Remote control drape

E

Irrigation sleeves x 2

F

Plug

G


Test chamber

H

Needle wrench

I

MICS angled needle



BL5612

20G High Speed (30-5000cpm)
Anterior Vitrectomy Cutter Pack

=6

STERILE



Phaco Tips for MICS® Surgery

RIS: Recommended Incision Size in Coaxial
OD: Outer Diameter
ID: Inner Diameter

Single Use

BL3318S MICS® needle 1,8   = 6



BL3318AS MICS® needle 1,8 angled 30°   = 6



RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm

Reusable limited re-use (10 cases)

BL3318 MICS® needle 1,8  = 6

BL3318A MICS® needle 1,8 angled 30°  = 6

RIS: 1.8-2.00mm OD: 0.74-0.95mm ID: 0.5-0.79mm

BL3118 Kit of 2 SU MICS Sleeves   = 12



For Standard Surgery

RIS: Recommended Incision Size in Coaxial
OD: Outer Diameter
ID: Inner Diameter

Single Use

DP8130S Standard needle 30° 19G   = 6

DP8145S Standard needle 45° 19G   = 6



RIS: 3.00-3.20mm OD: 1.70mm ID: 0.90mm

DP8230S Microflow needle 30°   = 6

DP8245S Microflow needle 45°   = 6



RIS: 2.50-2.65mm OD: 1.07mm ID: 0.90-0.51mm

DP8345S Microflow Plus needle 45°   = 6

DP8330S Microflow Plus needle 30°   = 6



RIS: 2.75-2.85mm OD: 1.17mm ID: 0.98-0.61mm

DP8430S Thin Tip needle 30° 20G   = 6



RIS: 2.5-2.65mm OD: 0.90mm ID: 0.67mm

DP8730S Standard needle 10°   = 6



RIS: 2.2mm OD: 0.90mm ID: 0.70-0.50mm

Reusable limited re-use (10 cases)

DP8130 Standard needle 30° 19G  = 6

DP8145 Standard needle 45° 19G  = 6



RIS: 3.00-3.20mm OD: 1.07mm ID: 0.90mm

DP8200 Microflow needle 0°  = 6

DP8215 Microflow needle 15°  = 6

DP8230 Microflow needle 30°  = 6

DP8230A Microflow needle 30° angled  = 6

DP8245 Microflow needle 45°  = 6



RIS: 2.50-2.65mm OD: 1.07mm ID: 0.90-0.51mm

DP8300 Microflow Plus needle 0°  = 6

DP8315 Microflow Plus needle 15°  = 6

DP8330 Microflow Plus needle 30°  = 6

DP8330A Microflow Plus needle 30° Angled  = 6

DP8345 Microflow Plus needle 45°  = 6



RIS: 2.75-2.85mm OD: 1.17mm ID: 0.98-0.61mm

DP8415 Thin Tip needle 15° 20G  = 6

DP8430 Thin Tip needle 30° 20G  = 6



RIS: 2.5-2.65mm OD: 0.90mm ID: 0.67mm

DP5531 kit of 2 SU Sleeves 2,2 surgery   = 12



DP5501 Standard kit of 2 SU Sleeves   = 12



Manipulators

SUF01 Koch Stop & Chop Manipulator   = 12



SUF02 Sinskey Hook   = 12



SUF03 Phaco Spatula   = 12



SUF04 Drysdale Nucleus Manipulator   = 12



SUF05 Horizontal Manipulator with Y-Hook   = 12



SUF06 Vertical Manipulator with Y-Hook   = 12



Cannulas

82020S Corneal Irrigator   = 60



Speculum

SUH01 Barraquer Eye Speculum   = 24



SUH04 Eye Speculum, adjustable   = 6



Forceps

SUE01-12 Capsulorhexis Forceps   = 12



SUE02-12 Bonn Forceps   = 12



SUE03-12 Tying Forceps   = 12



SUE04-12 McPherson Tying Forceps, angled   = 12



SUE05-12 Pierse Tying Forceps   = 12



SUE06-12 MICS™ Capsulorhexis Forceps, fine angled tips, angled shafts for 1.8mm incision.   = 12



SUE07-12 MICS™ Capsulorhexis with markings   = 12



SUE08-12 MICS™ Capsulorhexis with curved jaws   = 12



SUI20-12 MICS™ 1.0mm Capsulorhexis Forceps   = 12



Laseredge®



Clear Corneal Knives

Laseredge®
E7594 1.3mm Angled Double Bevel
A: 12.6mm B: 1.3mm C: 1.4mm




E7595 1.6mm Angled Double Bevel
A: 12.6mm B: 1.6mm C: 1.7mm

E7527 2.75mm Angled Double Bevel
A: 14.4mm B: 2.7mm C: 3.1mm

E7528 2.85mm Angled Double Bevel
A: 14.4mm B: 2.8mm C: 3.1mm

E7529 3.0mm Angled Double Bevel
A: 14.4mm B: 3.0mm C: 3.2mm



Slit Knives

Laseredge®
E7548ADB 2.2mm Angled Double Bevel
A: 14.4mm B: 2.2mm C: 2.7mm



E7549ADB 2.4mm Angled Double Bevel
A: 14.4mm B: 2.4mm C: 2.9mm

E7550A 2.5mm Angled Bevel Up
A: 14.4mm B: 2.5mm C: 3.8mm

E7550ADB 2.6mm Angled Double Bevel
A: 14.4mm B: 2.6mm C: 3.3mm

E7551A 2.65mm Angled Bevel Up
A: 14.6mm B: 2.6mm C: 3.9mm


E7559A 2.75mm Angled Bevel Up
A: 14.4mm B: 2.7mm C: 3.2mm

E7552ADB 2.80mm Angled Double Bevel
A: 14.4mm B: 2.8mm C: 3.4mm

E7553A 2.85mm Angled Bevel Up
A: 14.1mm B: 2.8mm C: 3.3mm

E7556A 3.0mm Angled Bevel Up
A: 14.4mm B: 3.0mm C: 3.9mm

E7557A 3.2mm Angled Bevel Up
A: 14.4mm B: 3.2mm C: 3.9mm



Stab Knives


Laseredge®
E7515 15°
A: 11.2mm B: 1.6mm C: 5.9mm



E7522 22.5°
A: 10.7mm B: 1.7mm C: 4.0mm

E7530 30°
A: 10.4mm B: 1.7mm C: 2.9mm

E7545 45°
A: 8.9mm B: 1.9mm C: 1.9mm



Trapezoidal Knives

Laseredge®
E7599 1.5-1.7mm Angled
A: 14.4mm B: 1.7mm C: 2.0mm D: 1.5 E: 2.9




E7596 1.5-2.0mm Angled
A: 14.4mm B: 2.0mm C: 1.8mm D: 1.5 E: 3.4

E7600 1.6-1.8mm Angled
A: 14.4mm B: 1.8mm C: 1.8mm D: 1.6 E: 2.9

E7602 1.8-2.2-2.8mm Angled
A: 16.8mm B: 2.2mm C: 1.5mm D: 1.8 E: 3.4

E7603 1.85-2.2mm Angled
A: 14.4mm B: 2.2mm C: 1.8mm D: 2.0 E: 3.7

E7601 2.0-2.2mm Angled
A: 14.4mm B: 2.2mm C: 2.3mm D: 2.0 E: 3.1



Stellaris® System Coaxial I/A Handpieces

85794ST MICS® I/A Handpiece, curved
shaft with sleeve for 1.8mm incision



85795ST MICS® I/A Handpiece, 45°
angled tip with sleeve for 1.8mm incision



85783ST 17G (Ø 1.5mm) I/A Handpiece,
curved shaft



85784ST 17G (Ø 1.5mm) I/A Handpiece
straight shaft 45° angled tip, 0.3mm
aspiration port



85786ST 16G (Ø 1.65mm) I/A
Handpiece straight shaft, 90° angled tip,
0.3mm aspiration port



Stellaris® System CapsuleGuard® I/A Coaxial
Handpieces

85910ST CapsuleGuard® I/A Handpiece
45° Tip – Incisions: 2.2mm-2.8mm



85912ST CapsuleGuard® I/A Handpiece
Straight Tip – Incisions: 2.2mm-2.8mm



85913ST CapsuleGuard® I/A Handpiece
Curved Tip – Incisions: 2.2mm-2.8mm



85914ST CapsuleGuard® I/A Handpiece
MICS® Straight Tip – Incisions: sub 20mm



85915ST CapsuleGuard® I/A Handpiece
MICS® 45° Tip – Incisions: sub 20mm



Stellaris® System Biaxial I/A Handpieces
Aspiration

- 85901S19 Aspiration Handpiece 19G
(Ø 1.07mm) STERILE  = 12
- 85901S20 Aspiration Handpiece 20G
(Ø 0.89mm) STERILE  = 12
- 85780S Aspiration Handpiece 21G
(Ø 0.8mm) STERILE  = 12
- 85790S Aspiration Handpiece 22G
(Ø 0.7mm) STERILE  = 12
- 85792S Aspiration Handpiece 23G
(Ø 0.6mm) STERILE  = 12



curved shaft roughened tip

Stellaris® System Biaxial I/A Handpieces Irrigation

- 85902ST19 1.8mm B-MICS 19G
(Ø 1.07mm) Irrigation Handpiece STERILE  = 12
- 85787ST 21G Irrigation Handpiece
(Ø 0.8mm) STERILE  = 12
- 85806ST 22G Irrigation Handpiece
(Ø 0.7mm) STERILE  = 12
- 85807ST 23G Irrigation Handpiece
(Ø 0.6mm) STERILE  = 12



curved shaft smooth tip

Stellaris® System Irrigating Chopper for biaxial
MICS

- 85903ST19 19G (Ø 1.07mm) Nagahara
type tip, end out STERILE  = 12



Stellaris® System Coaxial Zero Phaco
Handpieces

- SUP01ST Stellaris® System Zero Phaco
Handpiece 1.8mm STERILE  = 12
- For aspirating softened nuclei during femtosecond cataract surgery. Designed with an open port aspiration tube with 30° bevel. Recommended for use in a 1.8 mm incision. Designed for use with the Stellaris® System.



- SUP06ST Stellaris® System Zero Phaco
Handpiece 2.2mm STERILE  = 12
- For aspirating softened nuclei during femtosecond cataract surgery. Designed with an open port aspiration tube with 30° bevel. Recommended for use in a 2.2mm incision. Designed for use with the Stellaris® System.



- SUP05 Zero Phaco Aspiration
Handpiece 21G (0.8mm). STERILE  = 12
- Bimanual aspiration handpiece, designed with an open port aspiration tube for aspirating softened nuclei during Femtosecond cataract surgery. 30° tapered tip. For use with 85787ST - Irrigation Handpiece 21G.



Product Reference	Product Name	Page
18060	I-Pack™	24
18061	I-Pack™	24
18064	I-Pack™	24
18067	I-Pack™	24
18069	I-Pack™	24
09.05.23PIN	Pick Forceps 23G	23
09.05.25PIN	Pick Forceps 25G	23
09.10.23PIN	Fine Tip Eckardt Forceps 23G	22
09.10.23PIN-M	Fine Tip Eckardt Forceps 23G, 37mm	23
09.10.25PIN	Fine Tip Eckardt Forceps 25G	22
09.10.25PIN-M	Fine Tip Eckardt Forceps 25G, 37mm	23
09.12.23PIN	End-Grasping Forceps 23G	22
09.12.25PIN	End-Grasping Forceps 25G	22
09.12PIN	End-Grasping Forceps 20G	22
10.24.23PIN	Morris Vertical Scissors 23G	23
10.24.25PIN	Morris Vertical Scissors 25G	23
10.24PIN	Morris Vertical Scissors 20G	23
11.04.23PIN	Pyramid Tip Forceps 23G	22
11.04.25PIN	Pyramid Tip Forceps 25G	22
11.06.23PIN	Serrated Jaw Forceps 23G	22
11.06.25PIN	Serrated Jaw Forceps 25G	22
11.20PIN	Foreign Body Forceps 19G	23
11.40PIN	BB Removal Forceps 18G	23
12.06.25	Hubbard Visco Dissector 25G, 28G Tip	30
12.06.27	Hubbard Visco Dissector 27G	30
12.21PIN	Foreign Body Magnet 18G	23
20.04	Tano 20G	30
20.04.23	Tano 23G	30
20.04.25S (Stiff)	Tano 25G	30
20.04.27	Tano 27G	30
20.07	Stiff Tano 20G	30
20.16	Retractable Tano 20G	30
20.16.23	Retractable Tano 23G	30
20.16.25	Retractable Tano 25G	30
20.18	Adjustable Tano with Pick 20G	30
23.76.23PIN	Horizontal Scissors Curved 23G	23
23.76.25PIN	Horizontal Scissors Curved 25G	23
23.76PIN	Horizontal Scissors Curved 20G	23
30.32	6mm Tip Disposable Infusion Cannula 20G	30
30.32.23	6mm Tip Disposable Infusion Cannula 23G	30
30.32.25	6mm Tip Disposable Infusion Cannula 25G	30
30.90	Silicone oil collector chamber	31
39.03.23PIN	Tano Asymmetrical Micro Forceps 23G	22
39.03.25PIN	Tano Asymmetrical Micro Forceps 25G	22
39.03.27PIN	Tano Asymmetrical Micro Forceps 27G	22
39.07.23PIN	Apte Corset Forceps 23G	22
39.07.25PIN	Apte Corset Forceps 25G	22
39.08.23PIN	Eckardt Forceps 23G	22
39.08.25PIN	Eckardt Forceps 25G	22
39.08.27PIN	Eckardt Forceps 27G	22
39.08PIN	Eckardt Forceps 20G	22
39.12.23PIN	Microserrated™ End-Grasping Forceps 23G	22
39.12.25PIN	Microserrated™ End-Grasping Forceps 25G	22
39.15.23PIN	Microserrated™ Tewari Omnicap Forceps 23G	22
39.15.23PIN-M	Microserrated™ Tewari Omnicap Forceps 23G, 37mm	23
39.15.25PIN	Microserrated™ Tewari Omnicap Forceps 25G	22
39.15.27PIN	Microserrated™ Tewari Omnicap Forceps 27G	22
39.15.25PIN-M	Microserrated™ Tewari Omnicap Forceps 25G, 37mm	23
40.02-2	Disposable Iris Retractor	30

Product Reference	Product Name	Page
40.02-6	Disposable Iris Retractor	30
41.00	Awh Sutureless Lens Ring Disposable Lens	25
53.00.23	Bipolar Cautery 23G	30
53.00.25	Bipolar Cautery 25G	30
53.00.27	Bipolar Cautery 27G	30
55.12	Curved Laser Probe, 22° Curve 20G	16
55.21	Straight Laser Probe 20G	16
55.21.23	Straight Laser Probe 23G	16
55.21.25S	Straight Laser Probe 25G	16
55.21.27	Straight Laser Probe 27G	16
55.22	Straight Maxillum™ Illuminated Laser Probe 20G	17
55.23	TruCurve Curved Laser Probe, 37° Curve 20G	16
55.23.23	TruCurve Curved Laser Probe, 37° Curve 23G	16
55.23.25	TruCurve Curved Laser Probe, 37° Curve 24G	16
55.23.27	TruCurve Curved Laser Probe, 37° Curve 25G	16
55.25.23	Flexible Tapered Laser Probe 23G	16
55.25.25	Flexible Tapered Laser Probe 25G	16
55.26	Directional™ Laser Probe 20G	16
55.26.23	Directional™ Laser Probe 23G	16
55.26.231	Laser Probe Directional™ I 23G	16
55.26.23E	Extendable Directional™ Laser Probe 23G	16
55.26.25	Directional™ Laser Probe 25G	16
55.26.251	Laser Probe Directional™ I 25G	16
55.26.25E	Extendable Directional™ Laser Probe 25G	16
55.26.25S	Directional™ Laser Probe 25G	16
55.26E	Extendable Directional™ Laser Probe 20G	16
55.26.27	Directional™ Laser Probe 27G	16
55.27	Directional™ II Laser Probe 20G	16
55.27.23	Directional™ II Laser Probe 23G	16
55.27.25	Directional™ II Laser Probe 25G	16
55.27.27	Directional™ II Laser Probe 27G	16
55.29	Curved Maxillum™ Laser Probe, 37° Curve 20G	17
55.36.23E	Inverted Directional™ Laser Probe 23G	16
55.36.25E	Inverted Directional™ Laser Probe 25G	16
55.47.23P	Illuminated Directional™ Laser Probe 23G	17
55.47.25P	Illuminated Directional™ Laser Probe 25G	17
55.47P	Illuminated Directional™ Laser Probe 20G	17
55.48.23P	Illuminated Directional™ II Laser Probe 23G	17
55.48.25P	Illuminated Directional™ II Laser Probe 25G	17
55.62.23P	Straight Fixed Extended Illuminated Laser Probe 23G	17
55.62.25P	Straight Fixed Extended Illuminated Laser Probe 25G	17
55.69.23P	Curved Fixed Extended Illuminated Laser Probe, 37° Curve 23G	17
55.69.25P	Curved Fixed Extended Illuminated Laser Probe, 37° Curve 25G	17
55.69P	Curved Fixed Extended Illuminated Laser Probe, 37° Curve 20G	17
55.70.23P	Flexible Tapered Illuminated Laser Probe 23G	17
55.70.25P	Flexible Tapered Illuminated Laser Probe 25G	17
56.02.23P	Focal 23G	13
56.02.25PS	Focal 25G	13
56.02.27P	Focal 27G	13
56.02P	Focal 20G	13
56.07.23P	Chow Membrane Pick 23G	13
56.12.23P	Corona Shielded Widefield 23G	13
56.12.25P	Corona Shielded Widefield 25G	13
56.14	High-Flow Irrigating Endoilluminator 19G	13
56.20.23P	Vivid 23G	13
56.20.25P	Vivid 25G	13
56.21.23P	Widefield Diffusion 23G	12
56.21.25PS	Widefield Diffusion 25G	12
56.21P	Widefield Diffusion 20G	12

Product Reference	Product Name	Page
56.22.25P	Midfield 25G	13
56.22.27P	Midfield 27G	13
56.22P	Midfield 20G	13
56.24.23P	Aspirating Endoilluminator 23G	13
56.30.23P	Standard Illuminated Infusion Chandelier 23G	13
56.30.25P	Standard Illuminated Infusion Chandelier 25G	13
56.30P	Standard Illuminated Infusion Chandelier 20G	13
56.34.25P	High-Flow Illuminated Infusion Chandelier 25G	13
56.50.25P	Chandelier Awh 25G	13
56.50.29P	Chandelier Oshima Dual 29G	13
56.54.25P	Chandelier Awh Vivid 25G	13
56.54.27P	Chandelier Oshima Vivid 27G	13
57.20	Synergetics Illumination Adapter / Alcon Accurus, Xenon	12
57.26	Synergetics Illumination Adapter / DORC Associate	12
57.27	Synergetics Illumination Adapter / DORC Brightstar	12
57.30	Synergetics Illumination Adapter / Alcon Constellation	12
57.40	Synergetics Illumination Adapter / B+L Stellaris® PC	12
58.20	ONE-STEP™ 20G Polymide Cannula	11
58.20V	Valved ONE-STEP™ Port Entry System 20G	11
58.23	ONE-STEP™ 23G Polymide Cannula	11
58.23S	ONE-STEP™ 23G	11
58.23STH	Stealth Port Entry System 23G	11
58.23STHV	Valved Stealth Port Entry System 23G	11
58.23SV	Valved ONE-STEP™ 23G	11
58.25S	ONE-STEP™ 25G	11
58.25STH	Stealth Port Entry System 25G	11
58.25STHV	Valved Stealth Port Entry System 25G	11
58.25SV	Valved ONE-STEP™ 25G	11
58.27S	ONE-STEP™ 27G	11
58.27SV	Valved ONE-STEP™ 27G	11
65023S-12	Endgripping Forceps 23G	19
65025S-12	Endgripping Forceps 25G	19
82020S	Corneal Irrigator	35
85780S	Aspiration Handpiece 21G	38
85783ST	17G I/A Handpiece Curved	37
85784ST	17G I/A Handpiece Straight	37
85786ST	16G I/A Handpiece Straight	37
85787ST	21G Irrigation Handpiece (Ø 0.8mm)	38
85790S	Aspiration Handpiece 22G	38
85792S	Aspiration Handpiece 23G	38
85794ST	MICS® I/A Handpiece, Curved Shaft	37
85795ST	MICS® I/A Handpiece, 45° Angled Tip	37
85806ST	22G Irrigation Handpiece (Ø 0.7mm)	38
85807ST	23G Irrigation Handpiece (Ø 0.6mm)	38
85901S19	Aspiration Handpiece 19G	38
85901S20	Aspiration Handpiece 20G	38
85902ST19	1.8mm B-MICS 19G (Ø 1.07mm) Irrigation Handpiece	38
85903ST19	19G (Ø 1.07mm)	38
85910ST	CapsuleGuard® I/A Handpiece 45° Tip	37
85912ST	CapsuleGuard® I/A Handpiece Straight Tip	37
85913ST	CapsuleGuard® I/A Handpiece Curved Tip	37
85914ST	CapsuleGuard® I/A Handpiece MICS® Straight Tip	5
85915ST	CapsuleGuard® I/A Handpiece MICS® 45° Tip	37
BL14304	Stellaris® PC Next Generation	4
BL14334	Stellaris® PC Next Generation	4
BL2272	Laser Upgrade Kit	5
BL2272	Laser Upgrade Kit	14
BL2295	Integrated Foot Pedal	4
BL2295	Integrated New Design Footpedal	14

Product Reference	Product Name	Page
BL2296	Auxiliary Footpedal for LIO Use	15
BL2298	Dual LIO	15
BL3118	Kit of 2 SU MICS Sleeves	34
BL3170	Stellaris Phaco Handpiece	5
BL3234	Laser On/Off Keys (1 set)	14
BL3235	Laser Safety Goggles	14
BL3236	Zeiss Fixed ESF	14
BL3237	Leica Fixed ESF	14
BL3239	Laser Filter ESF Mount Topcon	14
BL3240	Zeiss 2 Position Filter ESF	14
BL3241	Leica 2 Position Filter ESF	14
BL3242	Adaptor for Fitting of 2 Position Filters	14
BL3244	Interlock Bypass Connector Kit	14
BL3245	Moeller Weidel Fixed ESF	14
BL3246	LIO Replacement Bulb	14
BL3247	Leica 4-way Filter	14
BL3270	Fragmentation Handpiece	5
BL3318	MICS® Needle 1,8 Reusable	34
BL3318A	MICS® Needle 1,8 Angled 30° Reusable	34
BL3318AS	MICS® Needle 1,8 Angled 30° Single Use	34
BL3318S	MICS® Needle 1,8 Single Use	34
BL3320	Xenon Lamp Replacement	5
BL3321	Mercury Lamp Replacement	5
BL3394	Handle Extension	10
BL4363	Stellaris® PC Zero Level Bottle Hanger - Dual Bottle	4
BL4364	Stellaris® PC/Stellaris Mayo Tray Adapter	5
BL4390	Foot Control Battery	4
BL4391	Foot Control Wall Charger	4
BL4393	Foot Control Charging Cradle	4
BL4394	Foot Control Power Cable	4
BL5110	Stellaris® Premium Vacuum Phacoemulsification Pack	32
BL5111	Stellaris® Basic Vacuum Phacoemulsification Pack	32
BL5112	Stellaris® StableChamber™ Vacuum Pack	32
BL5113	Stellaris® StableChamber™ Phaco Pack. MICS® Straight Needle	32
BL5114	Stellaris® Premium Vacuum Phaco Pack. MICS® Straight Needle	33
BL5114A	Stellaris® Premium Vacuum Phaco Pack. MICS® Angled Needle	33
BL5234	Adjustable Chandelier 23G	12
BL5235	Adjustable Chandelier 25G	12
BL5240	27G ESA Kit-Valved	10
BL5241	27G Single Trocar-Valved	10
BL5275	20G Fragmentation Needle 15°	5
BL5277	23G Fragmentation Needle 30°	5
BL5281	23G ESA Kit-Valved	10
BL5282	25G ESA Kit-Valved	10
BL5283	23G Single Trocar-Valved	10
BL5284	23G Replacement Valves	10
BL5285	25G Single Trocar-Valved	10
BL5286	25G Replacement Valves	10
BL5290	20G STD. Laser PROBE Straight	15
BL5290ASP	20G Soft Tip Asp. Laser PROBE S	15
BL5290ASPH	20G Hard Tip Asp. Laser PROBE	15
BL5290C	20G STD. Laser PROBE Curved	15
BL5290L	20G ILL. Laser PROBE Straight	15
BL5290LC	20G ILL. Laser PROBE Curved	15
BL5293	23G STD. Laser PROBE Straight	15
BL5293ASP	23G Soft Tip Asp. Laser PROBE S	15
BL5293ASPH	23G Hard Tip Asp. Laser PROBE	15
BL5293C	23G STD. Laser PROBE Curved	15
BL5293L	23G ILL. Laser PROBE Straight	15

Product Reference	Product Name	Page
BL5293LC	23G ILL. Laser PROBE Curved	15
BL5295	25G STD. Laser PROBE Straight	15
BL5295C	25G STD. Laser PROBE Curved	15
BL5295L	25G ILL. Laser PROBE Straight	15
BL5295LC	25G ILL. Laser PROBE Curved	15
BL5400	Basic Posterior Pack for procedural customization	9
BL5420W	20G Posterior Standard Pack	9
BL5423WV	23G Posterior Standard Pack	9
BL5425WV	25G Posterior Standard Pack	9
BL5420WX	20G Posterior Pack	8
BL5423WVX	23G Posterior Pack	8
BL5425WVX	25G Posterior Pack	8
BL5500	Basic Combined Pack for procedural customization	7
BL5520W	20G Combined Standard Pack	7
BL5523WV	23G Combined Standard Pack	7
BL5525WV	25G Combined Standard Pack	7
BL5520WX	20G Combined Pack	6
BL5523WVX	23G Combined Pack	6
BL5525WVX	25G Combined Pack	6
BL5612	High Speed Anterior Vitrectomy Cutter Pack 20GA	33
BL5620	20G Vitrectomy Cutter	10
BL5623S	23G Vitrectomy Cutter	10
BL5625S	25G Vitrectomy Cutter	10
BL5627S	27G Vitrectomy Cutter Bi-Blade®	10
BL5720	Focal 20G	12
BL5723	Focal 23G	12
BL5725	Focal 25G	12
BL5820	Wide Angle 20G	12
BL5823	Wide Angle 23G	12
BL5823M	Mid-Field 23G	12
BL5823D	Diffuse Wide Angle 23G	12
BL5825	Wide Angle 25G	12
BL5825M	Mid-Field 25G	12
BL5825D	Diffuse Wide Angle 25G	12
BL5827M	Mid-Field 27G	12
BL7620	High Flow VFI Cannula 20G	31
BL7623	High Flow VFI Cannula 23G	31
BL7625	High Flow VFI Cannula 25G	31
CX9400	Reusable Bipolar Cord with 2 Pin Connector	19
CX9430	Reusable Bipolar Cord with Lemosa Connector	19
CX5710	Viscous Fluid Extraction Kits	31
CX5720	Viscous Fluid Injection Kits for Oxane Range	31
DP5501	Standard Kit of 2 SU Sleeves	34
DP5531	Kit of 2 SU Sleeves 2,2 Surgery	34
DP8130	Standard Needle 30° 19G Reusable, Limited re-use (10 cases)	34
DP8130S	Standard Needle 30° 19G	34
DP8145	Standard Needle 45° 19G Reusable, Limited re-use (10 cases)	34
DP8145S	Standard Needle 45° 19G	34
DP8200	Microflow Needle 0°	34
DP8215	Microflow Needle 15°	34
DP8230	Microflow Needle 30°	34
DP8230A	Microflow Needle 30° Angled	34
DP8230S	Microflow Needle 30°	34
DP8245	Microflow Needle 45°	34
DP8245S	Microflow Needle 45°	34
DP8300	Microflow Plus Needle 0°	34
DP8315	Microflow Plus Needle 15°	34
DP8330	Microflow Plus Needle 30°	34
DP8330A	Microflow Plus Needle 30° Angled	34

Product Reference	Product Name	Page
DP8330S	Microflow Plus Needle 30°	34
DP8345	Microflow Plus Needle 45°	34
DP8345S	Microflow Plus Needle 45°	34
DP8415	Thin Tip Needle 15° 20G	34
DP8430	Thin Tip Needle 30° 20G	34
DP8430S	Thin Tip Needle 30° 20G	34
DP8730S	Standard Needle 10°	34
E7515	Laseredge® Stab Knife 15°	36
E7522	Laseredge® Stab Knife 22.5°	36
E7527	Clear Corneal Knife/ 2.75mm Angled Double Bevel	36
E7528	Clear Corneal Knife/ 2.85mm Angled Double Bevel	36
E7529	Clear Corneal Knife/ 3.0mm Angled Double Bevel	36
E7530	Laseredge® Stab Knife 30°	36
E7545	Laseredge® Stab Knife 45°	36
E7548ADB	Laseredge® Slit Knife/ 2.2mm Angled Double Bevel	36
E7549ADB	Laseredge® Slit Knife/ 2.4mm Angled Double Bevel	36
E7550A	Laseredge® Slit Knife/ 2.5mm Angled Bevel Up	36
E7551A	Laseredge® Slit Knife/ 2.65mm Angled Bevel Up	36
E7552ADB	Laseredge® Slit Knife/ 2.80mm Angled Double Bevel	36
E7553A	Laseredge® Slit Knife/ 2.85mm Angled Bevel Up	36
E7556A	Laseredge® Slit Knife/ 3.0mm Angled Bevel Up	36
E7557A	Laseredge® Slit Knife/ 3.2mm Angled Bevel Up	36
E7559A	Laseredge® Slit Knife/ 2.75mm Angled Bevel Up	36
E7594	Clear Corneal Knife/ 1.3mm Angled Double Bevel	36
E7595	Clear Corneal Knife/ 1.6mm Angled Double Bevel	36
E7596	Laseredge® Trapezoidal Knife/ 1.5-2.0mm Angled	36
E7599	Laseredge® Trapezoidal Knife/ 1.5-1.7mm Angled	36
E7600	Laseredge® Trapezoidal Knife/ 1.6-1.8mm Angled	36
E7601	Laseredge® Trapezoidal Knife/ 2.0-2.2mm Angled	36
E7602	Laseredge® Trapezoidal Knife/ 1.8-2.2-2.8mm Angled	36
E7603	Laseredge® Trapezoidal Knife/ 1.85-2.2mm Angled	36
E7924	Bipolar Eraser 25G, Tapered, 2-pin Connector	19
E7925	Bipolar Eraser 25G, Blunt, 2-pin Connector	19
E7926	Bipolar Eraser 25G, Tapered, Lemosá Type Connector	19
E7927	Bipolar Eraser 25G, Blunt, Lemosá Type Connector	19
E7928	Scleral Bipolar Eraser	19
E7929	Bipolar Eraser 20G	19
E7930	Bipolar Eraser 23G	19
E7931	Bipolar Eraser 23G, 2-pin Connector	19
E7932	Bipolar Eraser, Tapered, Lemosá type connector 27G	19
E7933	Bipolar Eraser, Tapered, 2-pin Connector 27G	19
LP7519	LaseredgePlus® 19G MVR Blade	20
LP7520	LaseredgePlus® 20G MVR Blade	20
P2-PLAMP	Replacement Bulb PHOTON II	12
PHOTON-2X	Dual Port Light Source	12
PLAMP	Replacement Bulb PHOTON	12
SUE01-12	Capsulorhexis Forceps	35
SUE02-12	Bonn Forceps	35
SUE03-12	Tying Forceps	35
SUE04-12	McPherson Tying Forceps, Angled	35
SUE05-12	Pierse Tying Forceps	35
SUE06-12	MICS™ Capsulorhexis Forceps	35
SUE07-12	MICS™ Capsulorhexis with Markings	35
SUE08-12	MICS™ Capsulorhexis with Curved Jaws	35
SUF01	Koch Stop & Chop Manipulator	35
SUF02	Sinskey Hook	35
SUF03	Phaco Spatula	35
SUF04	Drysdale Nucleus Manipulator	35
SUF05	Horizontal Manipulator with Y-Hook	35

Product Reference	Product Name	Page
SUF06	Vertical Manipulator with Y-Hook	35
SUH01	Barraquer Eye Speculum	35
SUH04	Eye Speculum, Adjustable	35
SUI03-12	Asymmetric Peeling Forceps 23G	19
SUI03-12	Asymmetric Peeling Forceps 23G	19
SUI05-12	Asymmetric Peeling Forceps 25G	19
SUI07-12	ILM Forceps 23G	19
SUI08-12	ILM Forceps 25G	19
SUI09-12	ILM Forceps 20G	19
SUI10-12	Serrated Forceps 20G	19
SUI11-12	Serrated Forceps 23G	19
SUI12-12	Serrated Forceps 25G	19
SUI13-12	Endgripping Forceps 27G	19
SUI14-12	ILM Forceps 27G	19
SUJ02	Passive Aspiration Handpiece. 1.5mm Soft Tip Cannula 25G	20
SUJ03	Passive Aspiration Handpiece. 1.5mm Soft Tip Cannula 23G	20
SUJ06	Passive Aspiration Handpiece. Blunt Needle 23G	20
SUJ07	Passive Aspiration Handpiece. Blunt Needle 25G	20
SUJ08	Backflush/Extrusion Handpiece. Blunt Tip 20G	20
SUJ09	Backflush/Extrusion Handpiece. Soft Tip Brush 20G	20
SUJ10	Backflush/Extrusion Handpiece. Soft Tip 20G	20
SUJ11	Backflush/Extrusion Handpiece. Blunt Tip 23G	20
SUJ12	Backflush/Extrusion Handpiece. Soft Tip brush 23G	20
SUJ13	Backflush/Extrusion Handpiece. Soft Tip 23G	20
SUJ14	Backflush/Extrusion Handpiece. Blunt Tip 25G	20
SUJ17	Backflush/Extrusion Handpiece. Blunt Tip 27G	20
SUK01	Membrane Pick 23G	19
SUK02	Membrane Pick 25G	19
SUP01ST	Stellaris® System Zero Phaco Handpiece 1.8mm	38
SUP05	Zero Phaco Aspiration Handpiece 21G (0.8mm)	38
SUP06ST	Stellaris® System Zero Phaco Handpiece 2.2mm	38
V20LC/SV20LC	Binocular Lens 20D Large	26
V20LCD	Single Use 20D Lens	25
V28LC	Binocular Lens 28D	26
V28LCD	Single Use 28D Lens	25
V3MIR	3 Mirror Laser lens	28
V30LC	Binocular lens 30D Large	26
V30PD10	30° Prism	25
V78C	Double Aspheric – 78D	27
V90C	90D	27
VAC	Area Centralis	28
VBCD10	Bi-Concave Disposable Lens	25
VCRLVIT	Central Retinal	28
VCRLVITSSV/SVCRLVITSSV	Central Retinal SSV	26
VDGTLWF	Digital Wide Field	27
VDVVIT	DynaView 156	28
VFD10	Flat Disposable Lens	25
VFLATSSVD	Flat, SSV	25
VG3	3 Mirror Lens Gonio	28
VG4HM	4 Mirrors Lens Gonio	27
VG4HMHAN2	3 Mirror High Mag Gonio	28
VG4LNF	4 Mirror Lens Gonio	28
VHRC	HR Centralis	27
VHRWF	High Resolution Wide Field	28
VHRXVIT	HRX Vit	28
VHRXVITSSV	HRX Vit SSV	29
VMD10	Magnifying Disposable Lens	25
VMQVIT/ SVMQVIT	MiniQuad®	26
VMQVITACS/SVMQVITACS	MiniQuad® ACS® Peripheral Indirect Vitreoretinal Procedures	26

Product Reference	Product Name	Page
VMQVITSSV	MiniQuad SSV	28
VMQVITSSVACS	MiniQuad ACS SSV VIT	29
VMQXLVIT/ SVMQXLVIT	MiniQuad® XL	26
VMQXLVITSSV	MiniQuad XL SSV	28
VQFL/SVQFL	QuadrAspheric®	26
VRL100	DK-Line® 5ml	31
VRL110	DK-Line® 7ml	31
VRL20	Okta-Line® 5ml	31
VRL500	Oxane® 1300, 10ml Syringe	31
VRL600	Oxane® 5700, 10ml Syringe	31
VRL700	Oxane® HD, 10ml Syringe	31
VS66	Field of view: 80°/96°	27
VSCA	Sterilization Tray Small	29
VSFNC	Super Field® NC Wide Field	26
VSPXL	SuperPupil® XL	27
VSQUAD160	SuperQuad® 160	26
VSRS2	Suture Ring	29
VVITHAN-LG	VitreoLens Handle	27
VVITHAN-MXQL	VitreoLens Handle	27
VWFD10	Widefield Disposable Lens	25



www.bauschsurgical.eu



@BauschSurgical



Bausch + Lomb Surgical



Please contact your local representative for more information on Bausch + Lomb retina products.
© 2018 Bausch & Lomb Incorporated.

EMEA_SU_C_VR_18_003

BAUSCH + LOMB
See better. Live better.